

# ACTIVE

# ClearBow™



The Active ClearBow is an ideal appliance for anterior correction of slight relapse, crowding, rotations and closing small diastemas. With anterior tooth set up, lingual acrylic trimmed to the new ideal tooth placement, elastic secured around the wilson loop to the distal portion of the labial wire allows alignment to occur.

This appliance can also be used as a dual purpose device when prescribed solely for retention. Simply, squeeze the buccal tubes to secure the labial wire and remove the elastics.



Contact us now for pricing.

**1-800-267-7982**  
**lab@orthodontic.ca**



Proudly serving the dental community since 1972, Rideau Orthodontic is committed to superior patient care, unmatched customer service and impeccable quality at a fair value. Save time, save money and improve peace of mind with Rideau Digital.

