



ALTERNATIVE LIGHTWIRE FUNCTIONALS



**THE ALF APPLIANCE IS MOST IDEAL FOR
STABILIZING CRANIAL ALIGNMENT**



- Fixed appliance producing gentle forces for optimal cranial alignment.
- Lingual fabrication designs offer inconspicuous treatment.
- Non-acrylic designs provide added tongue space.
- The ALF appliance is commonly combined with splint therapy and/or fixed brackets.
- Bands on 1st molars for patients in the mixed dentition.
- Expedite treatment with Indirect Tripping.

Numerous appliance variations available.
Please allow 5 laboratory days for ALF fabrication.

Continuing Education DVD programs available.

Rideau Orthodontic Mfg. Ltd.

1st Laboratory in Canada to fabricate the ALF

69 Beckwith Street, North
Smiths Falls, Ontario CANADA

(800) 267-7982


Tel: (613) 283-6841

Fax: (613) 283-5386

www.orthodontic.ca

orthodontic@orthodontic.ca

SPLINT DESIGNS



- ✓ ALF appliances are fabricated with Elgiloy wire.
- ✓ 24k gold solder joints to prevent discoloration/erosion.
- ✓ Bondable undercuts to improve appliance retention.

Serving the dental profession since 1972.

Made in Canada

