


QUAD HELIX

FIXED/REMOVABLE MAXILLARY APPLIANCE



- ✓ Maintain or develop the maxillary arch
- ✓ Auxiliary arms to tip incisors labially
- ✓ Stimulate buccal/lingual crown torque
- ✓ Used in conjunction with Fixed Bracket Therapy
- ✓ Best combined with a Low-Friction Bracket System
- ✓ Stimulates a cranial release
- ✓ Finger activate prior to insertion



➤ The Quad Helix is fabricated from the same Elgiloy wire as the ALF appliance.

Rideau Orthodontic Mfg. Ltd.

69 Beckwith Street, North
Smiths Falls, Ontario CANADA

(800) 267-7982

Tel: (613) 283-6841

Fax: (613) 283-5386

www.orthodontic.ca
orthodontic@orthodontic.ca

Fabrication Options

Plain or Straightwire Bands

➤ All attachments are **LASER** welded

Please specify modifications on your prescription.

Visit www.orthodontic.ca under Laboratory/Fixed for further information.

**A UNIQUE “CASE EVALUATION”
SERVICE AVAILABLE**

Serving the dental profession since 1972.

Made in Canada

