


WILSON[®] RLA

FIXED/REMOVABLE MANDIBULAR ARCH DEVELOPMENT



- ✓ Removable Lingual Arch
- ✓ Preserve the sacred lower deciduous molar space
- ✓ Upright retrusive lower incisors
- ✓ Maintain mandibular arch development
- ✓ Increase arch length
- ✓ Provide 1st molar anchorage

The Wilson[®] RLA is fabricated from the same Elgiloy wire as the ALF appliance.

Fabrication Options

Plain or Straightwire Bands

WILSON[®] RLA APPLIANCE
ADJUSTMENTS AVAILABLE
ON OUR WEBSITE.

Rideau Orthodontic Mfg. Ltd.

69 Beckwith Street, North
Smiths Falls, Ontario CANADA

(800) 267-7982

Tel: (613) 283-6841

Fax: (613) 283-5386

www.orthodontic.ca
orthodontic@orthodontic.ca

Wilson[®] Loop Activation Caution

Bond ridges above the lingual arch on the lingual of the lower laterals (above the wire) to prevent the molars from migrating mesially as the reciprocal force uprights and labializes the lower incisors.

Please specify modifications on your prescription.

Visit www.orthodontic.ca under Laboratory/Fixed for further information.

**A UNIQUE "CASE EVALUATION"
SERVICE AVAILABLE**

Serving the dental profession since 1972.

Made in Canada

