


BIONATOR



E&OE 2010

MIXED DENTITION CLASS II CORRECTION

- ✓ Early/mixed dentition
- ✓ Correct Anterior-Posterior Relationship
- ✓ Arch Development
- ✓ Control vertical
- ✓ Correct a Dual Bite
- ✓ Strengthen Musculature
- ✓ Encourage Lip Seal

**Ideal Adult Nighttime Retainer
for Post-TMJD treatment**

Rideau Orthodontic Mfg. Ltd.

69 Beckwith Street, North
Smiths Falls, Ontario CANADA

(800) 267-7982

Tel: (613) 283-6841

Fax: (613) 283-5386

www.orthodontic.ca
orthodontic@orthodontic.ca

- Bionator I : A-P & OPEN Vertical
- Bionator II: A-P & CLOSE Vertical
- Bionator III: A-P & Maintain Vertical
- Bionator 13: A-P & Md Vertical Only
- Orthopedic Corrector:
As above with 3 Expansion Screws for incremental A-P correction along with transverse development.

Please specify modifications on your prescription.
Visit www.orthodontic.ca under Laboratory/Functionals
for further information.

**A UNIQUE "CASE EVALUATION"
SERVICE AVAILABLE**

Serving the dental profession since 1972.

Made in Canada

