

RICK-A-NATOR



ANTERIOR-POSTERIOR CORRECTION OR MAINTENANCE

- ✓ Fixed appliance
- ✓ Correct Half-Tooth Class II
- ✓ Maintain Anterior/Posterior Correction
- ✓ Use simultaneously with fixed bracket therapy
- ✓ Increase Vertical with Vertical Box Elastics
- ✓ Worn 24/7 including meals

Ideal Phase II Appliance
following Phase I
TWIN BLOCK treatment

Rideau Orthodontic Mfg. Ltd.

69 Beckwith Street, North
Smiths Falls, Ontario CANADA

(800) 267-7982

Tel: (613) 283-6841

Fax: (613) 283-5386

www.orthodontic.ca
orthodontic@orthodontic.ca

CAUTION

➤ **Not recommended for
anterior open bite cases.**

Please specify modifications on your prescription.

Visit www.orthodontic.ca under Laboratory/Functionals for further information.

**A UNIQUE "CASE EVALUATION"
SERVICE AVAILABLE**

Serving the dental profession since 1972.

Made in Canada

