

# POSITIONER/SIMPLIFIER RX

For shipping instructions, please inquire by calling 1-800-267-7982

**Rideau  
Orthodontic  
Mfg.  
Ltd.**

69 Beckwith Street, North  
Smiths Falls, Ontario K7A 2B6  
P.O. Box 636, Ogdensburg, NY  
13669-0636 USA  
(613) 283-6841  
Toll Free: (800) 267-7982  
Fax: (613) 283-5386

EMMETT GRIFFITHS  
Representative

PLEASE PRINT Account # \_\_\_\_\_ or Doctor's Name \_\_\_\_\_

Address \_\_\_\_\_ City \_\_\_\_\_ Province \_\_\_\_\_ Postal Code \_\_\_\_\_  
 New Account  New Address or Phone Number State \_\_\_\_\_ Zip \_\_\_\_\_

Patient's Name \_\_\_\_\_ Patient's Age \_\_\_\_\_  Phone Me Concerning This Case

Date Sent \_\_\_\_\_ Date Required \_\_\_\_\_

**Please Send:**  Shipping Boxes  Wrapping  Shipping Labels  Record Model Boxes  Retainer Containers

**Rx Pads Required:**  Standard  Cosmetic Twin Block  3D\* Fixed/Removable™  Positioner/Simplifier  Straightwire

## GENERAL INSTRUCTIONS

- Do not reset teeth - fabricate positioner as specified.
- Construct standard set-up and positioner as specified.
- Construct gnathological set-up and positioner as specified.
- Construct diagnostic set-up as prescribed.
- Fabricate athletic sports guard.
- Detachable helmet strap.

## SET-UP INSTRUCTIONS

- Duplicate our models and return both sets.
- Remove all brackets and bands.
- Do not remove upper first molar bands.
- Allow for lower lingual 3-3 retainer.
- Reset all teeth to best advantage.
- Reset only teeth circled.

R	8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8	L
	8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8	

## SPACE CLOSURE:

- Close all spaces
- Close spaces as feasible
- Leave space between \_\_\_\_\_

## ARCH WIDTH:

- Upper:  Maintain  Constrict  Expand  
 Lower:  Maintain  Constrict  Expand

## ANTERIOR ROOT TORQUE:

- Upper:  Maintain  Lingual \_\_\_\_\_°  Labial \_\_\_\_\_°  
 Lower:  Maintain  Lingual \_\_\_\_\_°  Labial \_\_\_\_\_°

## OCCUSAL PLANE:


- Flat  Curve of Spee  As is

## HORIZONTAL OVERJET:

- As is  Reduce  Increase

## VERTICAL OVERBITE:

- As is  Reduce \_\_\_\_\_ mm  Increase \_\_\_\_\_ mm




Signature **X** \_\_\_\_\_ Licence \_\_\_\_\_ Phone (\_\_\_\_) \_\_\_\_\_

## POSITIONER INSTRUCTIONS

- MATERIAL:**  
 Vanguard (transparent vinyl)  
 Impak (clear, heat softened)  
 Rubber  White  Black

## TRIM



End Positioner Distal to  6's  7's  8's

## ARTICULATION


- Average hinge axis
- Hinge axis tracing  Lateral headplate
- Adjustable articulator  
 Sam II  Whipmix  Denar

## OPTIONS

- Moulded air holes
- Clasps between  Ball End  Loop End



- Socket liners on teeth circled



## CANADIAN SIMPLIFIER INSTRUCTIONS SILICONE MATERIAL:

- Light
- Medium
- Heavy

## CANADIAN SIMPLIFIER

- Inner & Outer Bow
- Inner Bow
- Supply Extraoral with elastics, please specify.

**SPECIAL INSTRUCTIONS:** \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_